

PYR-SL-R-TI

PR-SL-R-TI

Pyramid: Sliding & Rotating (Titanium) Product Information

The **PYR-SL-R-TI** features a 12.7 mm (1/2") slide off center, with closed ends on the base and can be oriented for either AP or ML slide. The PYR-SL-R-TI also features a 360° rotating pyramid while offering a low profile of only 16.25mm. This component is manufactured of titanium and is compatible with the European 4-hole pattern. The rotation is locked into place by tightening the 4 adjustment set screws on the interfacing connection component, such as PDI's PR-SL. Rotation is positively controlled in 3 degree increments with the use of interfacing radial serrations. The slide is locked at the same time the four adjustment set screws of the interfacing component are tightened. If attaching to a PDI lock, torque the M6 bolts to 12 Nm (9ft lb). **US Pat No. 6,033,440 and 6,231,618.**

Suggested L-Code:

L5910 – BK
L5920 – AK



PYR-SL-R-TI

SPECIFICATIONS

Added Height:	16.25mm (0.64 in)
Weight:	150 g (5 oz)
Load Rating:	120 Kg (265 lbs)
	K1-K4 Activity Level

The **PR-SL-R-TI** offers sliding (12 mm from center) The PR-SL-R-TI also features a 360° rotating pyramid receiver while offering a low profile of only 30.5 mm. This component is manufactured of titanium and is compatible with the European 4-hole pattern. The rotation is locked into place by tightening the 4 adjustment set crews on the receiver. The rotation is positively controlled in 3 degree increments with the use of interfacing radial serrations. The slide function is locked at the same time the four adjustment set screws of the receiver are tightened. If attaching to a PDI lock, torque the M6 bolts to 12 Nm (9ft lb). **US Pat No. 6,458,163.**

360°Rotation

M8 Set Screws

Suggested L-Code:

L5910 – BK
L5920 – AK



Sliding

SPECIFICATIONS

Added Height:	30.5 mm (1.3 in)
Weight:	171g (6.10oz)
Load Rating:	120 kg (265 lb)
	K1-K4 Activity Level

Warranty Information: PDI warrants its components to purchasers of the prosthetic components only. No warranties are made to consumers. All components are warranted for 30 days from original purchase date. This warranty applies only to manufacturing defects and does not cover: (1) any damage due to misuse, altering of the component, improper service or accidents; (2) any secondary costs that might be incurred in the exchange of components; (3) wear; normal, or due to incorrect selection of components. During the warranty period, defective parts will be replaced or repaired at no charge. Missing bolts will be replaced at retail prices. This warranty is in lieu of all other warranties, expressed or implied, including any implied warranty of merchantability or fitness for particular purpose. The remedies provided above are exclusive; in no event shall Prosthetic Design, Inc. be liable for any special, indirect, incidental, or consequential damages.

Returns: A restocking fee totaling 2% of retail price will be charged to customer should the customer wish to return any PDI product for credit. The part must be unused, and in its original, unopened packaging. In addition, the return must be initiated within 30 days of its purchase in order to qualify for the credit. If returned after 30 days of purchase date a 15% restocking fee will be charged.



Prosthetic Design, Inc., 700 Harco Drive, Clayton, OH 45315 USA
www.prostheticdesign.com Telephone: 800-459-0177
International Telephone: 937-836-1464 Fax: 937-832-5361

